
Pulse Oximeter Manual de Usuario YM403

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
Fabricante:
SHENZHEN YIMI LIFE TECHNOLOGY Co., Ltd.
302, Building C, Youlitong Technology
Industrial Plant N°56 Qingsong Road, Laokeng
Community, Longtitan Street Pingshan District,
518118 Shenzhen, China

Shanghai International Holding
Corp. GmbH
Eiffelstrasse 80, 20537
Hamburg, Germany

Pulse Oximeter User Manual

YM403

1. Operation Instructions

1.1 The product Operation Instructions

1. Installing two AAA batteries into battery cassette in correct polarities and cover it.
2. Plug one of fingers into rubber hole of the Oximeter (it is best to plug the finger thoroughly) nail surface upward, then releasing the clamp

3. Press the switch button once on front panel.
4. Your finger do not tremble during the Oximeter is working Your bod is not recommended in moving status.
5. Read correspondent date from display screen.

Note - Please use the medical alcohol to clean the rubber touching the finger inside of Oximeter, and clean the test finger using medical alcohol before and after each test (The rubber inside of the Oximeter belongs medical rubber which has no toxin and no harmful to the skin of human being). When your finger is plugged into the Oximeter, your nail surface must be upward.
Disposable batteries should be unloaded when used for a long interval.

Warning

- * The maximum application time for a PULSE OXIMETER PROBE at a single site is 4 hours.
- * The misapplication of PULSE OXIMETER PROBE with excessive pressure for prolonged periods can induce pressure injury.
- * Ambient light (including photodynamic therapy), physical movement (USERS and imposed motion), diagnostic testing, low perfusion; electromagnetic interference, HF SURGICAL EQUIPMENT; dysfunctional haemoglobin, there will be windage for SpO2 value.

1.2 Button operating Instructions

1. Put into two AAA batteries according to the instructions, the device will turn on automatically and display interface 1, then put into finger for measuring, if there is finger for detection and without operation, it power off automatically in 8s.
2. When there is battery, but the Pulse Oximeter power off, press the button, it will be opened.
3. During the measurement (there is a measurement signal and figure), press the button shortly, the interface can turn from interface 1 to interface 6 circularly. Press the button longly during the measurement, it can turn to settings menu interface.
4. In the Settings menu interface, press button longly, it can choose 8 options circularly, they are Alarm switch (Alarm Sound on / off), the add and subtract of oxygen lower limit value (SPO2 Lo +/-), PR Sound switch (PR Hi +/- PR LO +/-), PR Sound switch (BeeP) and the withdraw from the menu (Exit). The selected options foreground color will become Blue Square.
5. In the Settings menu interface, press button shortly, we can operate the value add and subtract or switch of current options, they are: The SpO2 default alarm lower limit value is 90, the adjustable range is from 70 to 95 and pressing the button shortly can do add and subtract of step. The PR default alarm upper limit value is 110, the adjustable range is from 110 to 250 and pressing the button shortly can do add and subtract of Steps. The PR default alarm lower limit value is 50, the adjustable range is from 15 to 95 and pressing the button shortly can do add and subtract of 5 steps limit value 90, PR upper limit value: 110, PR lower limit value: 50, PR sound off when we replaced the battery.
7. In the Settings menu interface, choose to Exit option and press button shortly, it will quit the setting menu interface.

2. Alarm

The Pulse Oximeter uses the following alarm indicators: audible alarm and visual alarms.

Note-To correctly identify visual alarms, always observe the Pulse Oximeter within 1 meter of its position.

2.1 Low power alarm

when battery power is at lowest level, the battery capacity indicates empty in Display, remind users of replacement of battery.

Note - When power is lost for less than or equal to 30s, the ALARM SETTINGS prior to the power loss shall be restored automatically.

2.2 Audible alarm

The audible alarm is buzzer alarm; the auditory sound pressure range is 45dB-70dB

2.3 SpO2 and PR alarm

- * Default alarming Value of PR: High Limit : 110BPM; Lower Limit : 50BPM.
- * Default alarming Value of SpO2 : 90%.
- * When it is under operation conditions, if the measurement of the PR and SpO2 value exceed the setting alarm limit and the alarm switch is on, it has a sound alarm and its long buzzer alarm.
- * When it is under operation conditions, if the PR sound switch is on, the measurement PR will go with a PR sound.

To test the alarm according the following:

- a) Connect the Pulse Oximeter to simulator;
- b) the simulator setting; SpO2 90%, PR 200bpm;
- c) the Pulse Oximeter alarm setting; SpO2 95%, PR High Limit 110 bpm;
- d) to verify the Pulse Oximeter alarm.

3. SpO2 Specifications

The update rate for the SpO2 value and pulse rate is typically 1 second. Data averaging and other signal processing on the displayed and transmitted data values of SpO2 and pulse rate is not more than 20 seconds. Depending on the magnitude of difference between the alarm limit and the displayed value, the alarm signal generation delay may be from 1 second to 20 seconds. The maximum alarm condition delay is 4 seconds, the Alarm alarm signal delay is 20 seconds, the average alarm condition is 2 seconds, the average alarm signal is 10 seconds.

Because pulse oximeter measurements are statistically distributed, only approximately two-thirds of pulse oximeter equipment measurements can be expected to fall within the Arms value measured by a CO-oximeter.

The signal adequacy is indicated by waveform, and it is NORMALIZED.

Note-The PULSE OXIMETER EQUIPMENT is calibrated to display FUNCTIONAL OXYGEN SATURATION It is not necessary to have a SpO2 calibration when the Pulse Oximeter is in use.

The following table describes the SpO2 specifications

Table 3 SpO2 Specifications

Parameter	Specification
SpO2 measurement range	5%-99%
SpO2 accuracy	± 2% in the range of 70%-99% Other scope is not defined
SpO2 Alarm range	Consistent with the display range, the set step length is 1%
Pulse rate range	30 bpm-250 bpm
Pulse rate accuracy	±1 bpm
Pulse rate Alarm range	Consistent with the display range, the set step length is 5bpm
Data update cycle	0.25s-2s
SpO2 PR average	8s
Pulse rate sound	PR tone
Peak Wavelength range	500 nm-100nm
Maximum optical output power	150 mW
Pulse rate display	Numeric
Power dissipation	in normal measurement, less than 50 mA, in power off state, less than 100 uA;


- (1) Sensor accuracy was obtained by performing controlled hypoxia studies on healthy, non-smoking adult volunteers (according to ISO 80601-2-61). The SpO2 readings have been compared to CO-oximeter measurements on arterial blood samples. To represent the general population, data from at least 10 subjects (male and female) with a wide range of skin color was taken to validate SpO2 accuracy. Accuracy data for Masimo sensors is included in the sensor package insert.
- (2) Information about wavelength ranges can be useful for clinicians performing photodynamic therapy.

4. Power Specifications

The following table describes the power specifications

Parameter	Specification
Battery	DC 3V AAA (X2) The Pulse Oximeter powers off if the battery power is almost depleted.

Battery Installation

5. Environmental Specifications

The following table describes the environmental specifications

Parameter	Specification
Temperature	Operating: 0 C-40 C Storage: -20 C + 55 C for the device
Relative humidity	Operating: 15% 0-80% (noncondensing) Storage: 10% 0-93% (noncondensing)
Brometric pressure	Operating: 59 kPa-107.4 kPa Storage: 22 kPa-107.4 kPa

6. Repair and Maintenance

- 1.Regular inspection to make sure that no obvious damage existed to affect the safety and performance of device.
2. No flammable substance, overtop or lower temperature and humidity existed
3. When the device is dabbled or there is hydraulic set existed, stop operating
- 4.When lower power capacity light in red pls replace the battery right away.
5. Pls clean the surface before applying for detection.
6. Pls take out the battery when device is not used for a period of time.
7. Pls dispose the battery according to the local statute.
8. If there is dust or dirt on surface, 75% density of medical alcohol can be used to clean the surface. Pls use dry fabric with little alcohol to can be alcohol permeates into the device.
9. The transportation and storage conditions are :
Temperature : -20C-55°C ; Humidity : 10% -93%

7. Troubleshooting

Trouble	Possible Reason	Solution
The SpO2 and pulse Rate display instable	1.The finger is not places inside enough. 2.The finger is shaking or the users is moving.	1. Place the finger properly and try again. 2. Let the users keep calm.
The device can not turn on	1.The batteries are drained or almost drained. 2.The batteries are not inserted properly. 3. The device's malfunction.	1. Change batteries 2. Reinstall batteries 3. Please contact the local service center.
The indicator light is off suddenly	1. The device willpower off automatically when it gets no signal for 8 seconds. 2. The batteries are almost drained.	1. Normal 2. Change batteries
The SpO2 and pulse Rate alarm	The value has exceeded the high alarm limit or below the low the alarm limit.	Check the users physiological status, and alarm limit set up.

8. Accessory

Fingertip Pulse Oximeter:1pcs
User Manual :1pcs

Date of last revision: 24/10/2022 V2

Manufacturer:
SHENZHEN YIMI LIFE TECHNOLOGY Co., Ltd.
302, Building C, Youlitong Technology
Industrial Plant N°56 Qingsong Road, Laokeng
Community, Longtan Street Pingshan District,
518118 Shenzhen, China

Shanghai International Holding Corp. GmbH
Eiffelstrasse 80, 20537
Hamburg, Germany

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