# Kontaktloses Infrarot Thermometer

Model: UFR-106

---

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**Fabricante/Manufacturer:**  
SHENZHEN URION TECHNOLOGY Co., Ltd.  
Floor 4-6th of Building D, Jiale Science & Technology  
Industrial Zone, N°3, Shenzhen

**Shanghai International Holding Corp. GmbH**  
Eiffestrasse 80, 20537  
Hamburg, Germany

User manual in Deutsch,  
Français, Português, Italiano.

# Infrared Non-Contact Thermometer

Model: UFR-106

---

## User manual

Thanks for purchasing our product.

To ensure the correct use of this product, please read this manual carefully and keep this manual properly so that you can read it at any time.

### 1. Safety Instructions

**Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use. Should you give this device to another person, it is vital that you also pass on these instructions for use.**

- Please operate according to this manual to avoid damage to the device.
- Used to measure body temperature on the forehead and the surface temperature of an object.
- Please follow the doctor's instructions and do not make self-judgment and treatment based on the measurement results. The measurement results are for reference only.
- Must be operated by adults, and be placed out of the reach of infants.
- Please do not let children use the body alone, as this may cause injury.
- Do not strongly collide, drop, step on, or vibrate this product.
- Do not use this device to measure any objects whose temperature exceeds 100.0°C, otherwise the measurement cannot be performed correctly.
- Do not use it after entering or leaving an environment with a large temperature difference immediately to avoid measuring temperature deviation.
- Do not disassemble or modify the medical infrared thermometer by yourself, otherwise it will not be able to measure correctly and affect the product warranty.
- Clean the thermometer according to the requirement of this manual before and after each use to avoid possible bacterial infections.
- Avoid strong electromagnetic interference when using, such as near mobile phones, microwave ovens, etc.
- Do not disassemble or modify the product by yourself, otherwise it may cause product damage, inaccurate measurement and invalid warranty.
- Never immerse the thermometer into water or other liquids.

## 2. Explanation of symbols

	<b>IMPORTANT</b> Follow the instructions for use! Non-observance of these instructions can result in serious injury or damage to the device.
	<b>WARNING</b> These warning notes must be observed to prevent any injury to the user.
	<b>CAUTION</b> These notes must be observed to prevent any damage to the device.
	<b>NOTE</b> These notes give you useful additional information on the installation or operation.
	Device classification: type BF applied part
	Information about the manufacturer
	Complies with waste electrical and electronic equipment directive
IP22	Protection against harmful ingress of water or particulate matter

## 3. Indication for Use

The infrared Non-contact thermometer, Model IR-FM03 is intended for the intermittent measurement of forehead temperature of human body. The device is indicated for use by people of all ages in the home.

## 4. Features

- Easy to hold and operate
- Two measurement mode: human body temperature and object temperature mode
- LCD display with backlight, and there are three kinds of backlight display (Green, Yellow, Red)
- 32 sets of memories, and it is easy to recall
- Switching between °C or °F
- Resolution is 0.1°C (0.1°F)
- Low battery display
- Prompt tone setting

## 5. Product Main Structure


It is mainly composed of ABS plastic shell, copper head structure, PCB circuit board, IC, resistance, capacitance, infrared sensor, LCD, buzzer, and battery line.

### 5.1 Components Description

- ① Measure button
- ② Battery cover
- ③ Infrared sensor
- ④ Set button
- ⑤ Up button
- ⑥ Down button
- ⑦ Display screen

### 5.2 Display

## 6. Package Contents

The factory standard packaging of this device includes the following contents:

No.	Parts	Quantity
1	Infrared Non-contact thermometer	1
2	User manual	1
3	Quality Card	1

Note: Before use, please check whether the host and accessories are complete according to the packing list. If you find that the above physical objects are inconsistent with the information, please contact the distributor.

## 7. Software version



V1.0

## 8. Method of Application

### 8.1 Battery Installation and Replacement


#### A. Battery Installation

- The detachable part of this device is the battery, and its replacement method is to use DC.3V (2xAAA dry battery) for replacement, otherwise it will cause measurement error.
- Before installing the batteries, please confirm whether the specifications of the batteries are consistent with the manufacturer's specifications.
- Be sure to use the specified model parts, if necessary, please contact the seller or manufacturer.

What to do	Install the battery	Remove the battery
Picture tips		
Steps	<ul style="list-style-type: none"> <li>• Turn the infrared thermometer over and push the battery cover in the direction of the arrow.</li> <li>• Install 2xAAA (3.0V) batteries according to the battery polarity marked inside the battery box.</li> <li>• After loading the battery as shown in the figure, slide the battery cover in place.</li> </ul>	<ul style="list-style-type: none"> <li>• Push off the battery cover and remove the battery according to the arrow instruction</li> <li>• Close the battery cover as shown in the figure.</li> </ul>

#### B. Battery Replacement

When the battery symbol as below on the screen becomes hollow, it indicates that the power is insufficient and the batteries need to be replaced. In the shutdown state, remove the old batteries, and install two new AAA batteries according to section A above to replace the batteries.

#### Warning:

- Do not disassemble the battery.

- Please use quality-guaranteed batteries or those specified by this product. Otherwise, the leakage will corrode the casing and metal shrapnel, so that the equipment cannot work.
- Do not mix new and old batteries, and batteries of different brands and types.
- Please install the battery according to the suggested polarity, or will cause the device to display no power or insufficient power, so that the device cannot work or even be damaged.
- Please immediately replace the battery with insufficient power when the infrared thermometer shows a low voltage symbol to indicate that it is within the low voltage range.
- Please remove the battery if it will not be used for a long time (more than 3 months) to avoid battery leakage affecting the normal use of the device.
- If the battery is charged, or is improperly placed with positive and negative electrodes, short circuited, put into water, hit, punctured, etc., it may cause leakage, rupture and explosion.
- Dispose waste batteries according to the requirements of local environmental protection regulation.


## 8.2 Preparation and Inspection Before Measurement

- 1) This product needs to be used indoors with an ambient temperature of 59°F/104°F (15°C/40°C).
- 2) Place the infrared thermometer in the operational environment for more than 30 minutes to make the temperature of the infrared thermometer consistent with the operational environment.
- 3) Check whether the temperature sensing head of the thermometer is clean before each use, and make sure that there is no dirt and foreign matter blocking at the infrared sensor.

\* When there is a difference between the temperature of the operational environment and the temperature of the infrared thermometer, or when the infrared thermometer is used nearby an air conditioner, it will not measure correctly.


## 8.3 Temperature Unit Settings

Short press the "Measure" button to turn on the thermometer, and long press the "M" button for 3 seconds, when the screen displays: Fl, select "+" and "-" keys to switch between degrees Celsius °C and Fahrenheit °F.


## 8.4 Switch between body temperature mode and object temperature mode

The product has two modes: body mode and object mode. Short press the "Measure" button to turn on the thermometer, and short press the "M" button to switch between body mode and object mode.

## 8.5 Turn On/Off the Prompt Tone

In the power-on state, press and hold the "M" button for 3 seconds enter into the setting mode. At this time, short press the "M" button until the display shows "F4", indicating the prompt tone setting state has been entered. Short press "+" or "-" to turn on/off the prompt tone.

## 8.6 High Temperature Alarm Setting

In the power-on state, press and hold the "M" button for 3 seconds enter into the setting mode. At this time, short press the "M" button until the display shows "F2", indicating the high temperature alarm setting has been entered. The default for high temperature alarm limit is 38.0°C. short press the "+" or "-" button to plus or minus the alarm limit and then press "M" button to complete the setting.

## 8.7 Body Temperature Measurement

- Press the "measure" button to start the device. At this time, the LED displays a yellow backlight and displays all segments.
- Press mode button to make infrared thermometer enter body temperature measurement mode. An icon **Body** indicating that the device have entered the body temperature measurement mode.
- Align the sensing probe to the center of the forehead at a distance of 3-5 cm, and press the measure button to start measuring the temperature.
- When you hear a short beep, the body temperature has been measured and the measurement is complete and the measurement results can be read.
- Three measurements are recommended. If the three measurements are different, select the highest temperature value.

Note:

- Do not drink, eat or exercise before or during the measurement.
- After continuous measurement, please wait at least 2 minutes before turning on/off the infrared thermometer.
- Do not collide, fall, trample or vibrate the product strongly.

- Before aiming at the target, please remove water, dust or dirt.
- When an infrared thermometer is aimed at the forehead, please clean the scanning area and remove hair or sweat to ensure the probe is clean.
- When measuring human body temperature, the center of forehead should be aligned at a distance of 3-5 cm to avoid contact between the sensor head and human body.

## 8.8 Object Temperature Measurement



- Press the "measure" button to start the device. At this time, the LED displays a yellow backlight and displays all segments.
- Press mode key to make infrared thermometer enter object temperature measurement mode. An icon **Surface temp** indicating that the device enters the object mode.
- Aim the sensing probe at the object to be measured at a distance of 3-5 cm, then press the measure button to start measuring the temperature.
- When you hear a short beep, the object temperature has been measured and the measurement is complete and the measurement results can be read.
- Three measurements are recommended. If the three measurements are different, select the highest temperature value.

## 8.9 View and clear Memory Values

### A. Read Memory Values

The product will automatically record the data obtained in each test. Up to 32 measurement data can be recorded. You can view these data by following the steps below:

- In the standby state, press the "+" button, the device can display the last measured temperature.
- Press the "+" button again, the last data will be displayed, and so on, and you can view the stored measurement data.
- If without pressing any key, it will exit this mode within 10 seconds and automatically shut down.

The infrared thermometer can automatically save the most recent 32 measurement values. If the memory times are full, the earliest measurement result will be deleted.

## B. Memory Clearance

In the power-on state, you can clear all the current memory by holding down the "+" key for 3 seconds.

## 8.10 The Method to Calibrate Your Thermometer

The thermometer is initially calibrated at the time of manufacture. If this thermometer is used according to the user manual, it is regular recommended to do the re-calibration every year to ensure the functioning correct. The re-calibration should be done by a qualified technician. Please send your thermometer to the dealer for re-calibration every year. However, if any time you question the accuracy of temperature measurements, please contact customer service.

Notes

Calibration mode is used by qualified organizations and personnel to check the accuracy of the sensor. Using this mode by ordinary users will lead to wrong results.

## 9. Reference temperature

### A. Normal temperature range at different measurement parts

Normal Temperatures According to Measurement Parts

Measurement part	Normal Temp °C	Normal Temp °F
Rectal	36.6 to 38	97.8 to 100.4
Oral	35.5 to 37.5	95.9 to 99.5
Axillary	34.7 to 37.3	94.4 to 99.1
Ear	35.8 to 38	96.4 to 100.4

The temperatura of the human body varies throughout the day. It can also be influenced by numerous external factors: age, sex, type and thickness of skin.

### B. Normal body temperature range at different ages





Normal Temperatures According to Age

Age	Temp °C	Temp °F
0-2 years	36.4 to 38.0	97.5 to 100.4
3-10 years	36.1 to 37.8	97.0 to 100.0
11-65 years	35.9 to 37.6	96.6 to 99.7
> 65 years	35.8 to 37.5	96.4 to 99.5

Attention: ASTM laboratory accuracy requirements in the display of 37 to 39°C (98 to 102°F) for Skin IR thermometers is  $\pm 0.3^{\circ}\text{C}(\pm 0.5^{\circ}\text{F})$ , whereas for mercury in-glass and electronic thermometers, the requirement per ASTM Standards E667-86 and 1112-86 is  $\pm 0.1^{\circ}\text{C}(\pm 0.2^{\circ}\text{F})$ . Manufacturer will make available upon request specific instructions for tests to determine the laboratory error, clinical bias and clinical repeatability.

When you need to clinical accuracy characteristics and procedures, please contact with the local distributors who can obtain it from our company.

## 10. Troubleshooting of Abnormal Measurement Result

Display identification	Significance	Possible causes and methods to resolve
	Excessive measurement temperature	Temperature mode: the measured target temperature is higher than 43°C.
	Measuring temperature too low	Temperature mode: the measured target temperature is below 32°C.
	Battery power low prompt	Battery power is insufficient and "Lo" is displayed on the screen. Please replace the battery immediately.
	System failure	Return to the seller, or contact with the local distributors to consult handling method.

\* If the error message cannot be eliminated, please contact the manufacturer.

## 11. Cleaning, Storage and Disposal

### Cleaning

- Clean the thermometer according to the requirements of this manual after each use.
- Use a soft, dry cloth to wipe and clean the device. Be careful not to wipe with volatile oil or thinner.
- Do not allow liquids to enter the product. The recommended cleaning frequency is once a week.
- You can use the absorbent cotton dipped in 70% alcohol to gently wipe the case and sensor head (except the panel and nameplate). The recommended disinfection cycle is once a month.
- Do not clean the device with corrosive detergent.

### Storage

- Do not store or only use the device in damp, high temperature, falling, dust, sun and other environments.
- Do not disassemble or modify this medical infrared thermometer without permission. If you replace the original parts with parts not provided by the manufacturer, it may cause measurement errors.
- Do not use in the environment with electromagnetic wave interference (such as mobile phones or microwave ovens, amateur radio, AM and FM radio broadcasts and TV broadcasts).

- Please take out the disposable battery when not using this product for a long time, so as not to cause the dry battery to heat, leak, rupture, etc., which may damage the medical infrared thermometer.
- This medical infrared thermometer is a sensitive measuring instrument, it is recommended to calibrate it every year or after mechanical impact (such as falling).

### Disposal

When the machine reaches the end of its service life, there may be inaccurate measurement; discarded medical infrared thermometers and accessories may cause environmental pollution.

- Please treat and recycle the products according to the relevant requirements of national environmental protection, and the used batteries shall be treated according to the relevant environmental protection regulations of the city.
- Please dispose of used batteries in accordance with the relevant environmental protection regulations of the city.

## 12. Product Specification

Items	Specification
Product name	Infrared Non-Contact Thermometer
Product model	IR-FM03
Power Supply Voltage	d.c.3.0V (2pcs AAA batteries)
Measurement location	Forehead
Temperature Sensor	Infrared Sensor
Display Mode	Digital LCD
Proper measuring distance range	3 to 5 cm
Measurement Range	Body Mode: 32.0°C~42.9°C (89.6°F~109.4°F) Object Mode: 0°C~80°C (32°F~176°F)
Measurement Accuracy	±0.2 °C: 35.0 °C ~42.0°C ±0.3°C: 32.0 °C ~34.9°C, 42.1°C ~42.9°C  ±0.4°F: 95.0°F ~107.6°F ±0.5°F: 89.6°F ~94.8°F, 107.8°F ~109.4°F,
Repeatability	≤0.2 °C (0.4 °F)
Display Resolution	0.1°C/°F
Measurement Time	1 second
Automatic shut-down	Approx. 8 seconds after last measurement has been taken
Mode of operation	Continuous operation*
Memory	32 sets memories
Backlight	Red backlight: 38.5°C≤T≤43.0°C, fever Yellow backlight: 37.5°C≤T≤38.5°C, moderate fever Green backlight: T <37.5°C normal

Housing material	ABS Plastic
Patient-contact button material	ABS Plastic
Sensor type	
Operational Environment	Temperature: -20°C ~ 40°C Relative Humidity: 15%~95%
Transport/Storage Condition	Temperature: -20°C ~ 55°C Relative Humidity: 15%~95% Atmospheric Pressure: 70kPa~106kPa
Safety Classification	Internal power type, BF type application
EMC Classification	Group one, Class B
Weight	49g (without battery)
Dimensions	158mm×40mm×42mm (L×W×H)
Expected Service Life	5 years

\* Continuous Operation - means the device can operating in normal use for an unlimited period of time without the specified temperatura being exceeded.

### 13. EMC Declaration

<b>Guidance and Statement on Electromagnetic Launching</b>		
The infrared thermometer is intended for use in the electromagnetic enviroment specified below. The purchaser or user of the infrared thermometer should ensure that it is used in this electromagnetic environment:		
<b>Launch test</b>	<b>Compliance</b>	<b>Electromagnetic Environment-Guidelines</b>
Radio-frequency emission CISRP 11	Group 1	Infrared thermometers use RF energy only for its internal functions. Therefore, its RF emissions are very low and may not cause any interference to nearby electronic equipment.
Radio-frequency emission CISRP 11	Class B	Infrared thermometers are suitable for use facilities, including homes and public low-voltage power supply environment.
Harmonic emission IEC61000-3-2	Not applicable	
Voltage fluctuation/flicker emission IEC61000-3-3	Not applicable	

<b>Guidance and Statement on Electromagnetic Immunity</b>			
The infrared thermometer is expected to be used in the electromagnetic environment specified below. The purchaser or user of the infrared thermometer should ensure that it is used in this electromagnetic environment.			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Coincidence level</b>	<b>Electromagnetic Environment - Guidelines</b>
Electrostatic Discharge (ESD) IEC61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	The floor should be Wood, concrete or tile, and if the floor is covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient burst IEC61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable	Network power quality must be a typical commercial or hospital environment.
Surge IEC61000-4-5	± 1 kV differential mode ± 2 kV common mode	Not applicable	Network power quality must be a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines	0 % UT; 0,5 cycle g) At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles at 0° 0 % UT; 250/300 cycle	Not applicable	Network power quality must be a typical commercial or hospital environment. If the equipment needs to keep running continuously during the interruption of network power supply, we recommended using uninterruptible power supply UPS for power supply.
Power frequency (50Hz/60Hz) magnetic field IEC/EN 61000-4-8	30A/m	30A/m	The power frequency magnetic field must be the level of a typical commercial or hospital environment.
Note: UT refers to the AC grid voltage before test voltage is applied.			

## 14. After-sales service

- The infrared thermometer has a free warranty of at least one year from the date of sale. A guarantee card filled out and sealed by the Distributor and granted must be issued during the warranty.

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